

Invacare®  
Solara®  
Limited  
Wheelchair



# Invacare® Solara® Limited Tilt-in-Space Wheelchair



**Yes, you can.™**

The Solara® Limited pediatric tilt-in-space wheelchair provides consumers with one of the most prescribed weight shifting tilts when funding options are limited. Based on the revolutionary Solara, the Solara Limited pediatric wheelchair provides a durable, rigid base along with an appropriate number of the most important options required to meet the consumers' needs.



Solara Limited shown in red with 70° swingaway front rigging, flip-up composite footplates, 6" x 2" semi-pneumatic casters, 12" composite wheels with pneumatic tires, conventional adjustable height full length arms, adjustable angle rigidizing handle bar, and transport tie down brackets.

#### The Solara Limited Base Chair Includes

- Rigid, rear wheel weight shift tilt-in-space frame with 45° of posterior tilt and up to 10° of anterior tilt
- 200 lb. weight capacity
- Seat widths (12"-14")\*
- Seat depths (15"-18")
- Transport tie down option\*\*
- Cantilever, non-locking bolt on arms
- Fixed back height (20" or 24")
- Three frame finishes (Wet Black, Red or Midnight Blue)
- 12" composite wheels
- 70° and 90° center pivot front rigging
- Flip-up composite footplates
- Attendant operated wheel locks
- Permanent axles
- Anti-tippers
- Adjustable angle back
- 17" or 18.25" front seat to floor height (factory set at 17")

\*Growable to 18" wide with additional parts.

#### Other Available Options Include

- Recliner back
- Larger seat widths and depths (20", 21" and 22")
- Conventional and cantilever adjustable height arms
- Flush or drop mounted seat pans
- 6" x 2" semi-pneumatic casters
- Pneumatic tires
- Flat-free inserts
- Elevating or articulating legrests
- Adjustable angle Flip-up footplates
- Heel and toe loops
- Legrest accessories
- Invacare seating and positioning products
- Rigidizers and drop base
- Rigidizing handle bars
- Push button seat positioning strap



Front and rear transport tie down brackets are mounted where caregivers can reach.\*\*

Transport locking pins should be used to lock gear rack when chair is transported by bus or van.\*\*

#### Specifications

Frame style	Tilt
Seat widths	12"-14"
Seat depths	15"-18"
Back height	20" or 24"
Back angle adjustment	90° to 120° in 10° increments
Tubing diameter	1"
Tilt range	0° to 45°
Weight capacity	200 lb.
Weight of chair	43 lb.
Limited warranty	Lifetime on frames and crossbraces

\*\*The Solara® wheelchair was crash tested with an Invacare Personal Back™10 and Infinity® seat cushion and conforms with the requirements of ANSI/RESNA WC/19 vol. 1 – Section 19.5.3 (Frontal Impact Test). Invacare has not performed testing with any other seating system. See Solara sell sheet (form no. 04-036) for ANSI/RESNA WC/19 sec. 6.2 Presale Literature Notification. Invacare recommends that wheelchair users NOT be transported in vehicles of any kind while in wheelchairs. As of this date, the Department of Transportation has not approved any tie down systems for transportation of a user while in a wheelchair in a moving vehicle of any type.



Yes, you can.™

Invacare Corporation  
www.invacare.com

#### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

#### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

Invacare, Infinity and Solara are registered trademark and Yes, you can., the blue medallion design, Personal Back and Quick Set and Solara are trademarks of Invacare Corporation. Specifications are subject to change without notification. Accessories and seating products shown in the photographs are not necessarily included in the base price of wheelchair. Refer to current price list located at [www.invacare.com](http://www.invacare.com)

© 2004

Invacare Corporation  
Form No. 04-034



"Invacare says 'Yes, you can,'" and I couldn't have said it better myself."

*Arnold Palmer*

## Free Manuals Download Website

<http://myh66.com>

<http://usermanuals.us>

<http://www.somanuals.com>

<http://www.4manuals.cc>

<http://www.manual-lib.com>

<http://www.404manual.com>

<http://www.luxmanual.com>

<http://aubethermostatmanual.com>

Golf course search by state

<http://golfingnear.com>

Email search by domain

<http://emailbydomain.com>

Auto manuals search

<http://auto.somanuals.com>

TV manuals search

<http://tv.somanuals.com>